

## In-Gel Digest Protocol

### Stock Solution:

1. 25 mM  $\text{NH}_4\text{HCO}_3$  (100 mg/50 ml)
2. 25 mM  $\text{NH}_4\text{HCO}_3$  in 50% ACN
3. 50% ACN/5% formic acid (may substitute TFA or acetic acid)
4. 12.5 ng/ $\mu\text{L}$  trypsin in 25mM  $\text{NH}_4\text{HCO}_3$  (freshly prepared)
5. 10 mM (1.5 mg/mL) DTT in 25 mM  $\text{NH}_4\text{HCO}_3$  (freshly prepared)
6. 55 mM (10 mg/mL) iodoacetamide in 25 mM  $\text{NH}_4\text{HCO}_3$  (freshly prepared)

### Procedure:

1. Dice each gel slice into small pieces (1 mm<sup>2</sup>) and place into 0.65 mL siliconized tube.
2. Add ~100 $\mu\text{L}$  (or enough to cover) of 25mM  $\text{NH}_4\text{HCO}_3$ /50% ACN and vortex for 10 min.
3. Using gel loading pipet tip, remove the supernatant and discard.
4. Repeat steps 3 and 4 once or twice.
5. Speed Vac the gel pieces completely dry (~ 20 min).
6. Add 25  $\mu\text{L}$  (or enough to cover) 10 mM DTT in 25 mM  $\text{NH}_4\text{HCO}_3$  to dried gels. Vortex and spin briefly. Allow reaction to proceed at 56°C for 30 min.
7. Discard the supernatant, add 25  $\mu\text{L}$  55 mM iodoacetamide to the gel pieces. Vortex and spin briefly. Allow reaction to proceed in the dark for 45 min at room temperature.
8. Discard supernatant. Wash gels with ~100  $\mu\text{L}$   $\text{NH}_4\text{HCO}_3$ , vortex 10 min, and spin.
9. Discard supernatant. Dehydrate gels with ~100 $\mu\text{L}$  (or enough to cover) of 25 mM  $\text{NH}_4\text{HCO}_3$  in 50% ACN, vortex 5 min, spin. Repeat one time.
10. Dry the gel pieces completely with speed vacuum (~20 min). Proceed with trypsin digest.
11. Add trypsin solution to just barely cover the gel pieces. Estimate the gel volume and add about 3 $\times$  volume of trypsin solution. This volume will vary from sample to sample, but on average ~5-25  $\mu\text{L}$  is sufficient. Rehydrate the gel pieces on ice or at 4°C for 10 min. Spin.
12. Add 25mM  $\text{NH}_4\text{HCO}_3$  as needed to cover the gel pieces.
13. Spin briefly and incubate at 37°C for 4 hours – overnight.
14. Transfer the digest solution (or extraction) into a clean 0.65 mL siliconized tube (one per sample).
15. To the gel pieces, add 30  $\mu\text{L}$  (enough to cover) of 50% ACN/5% formic acid, vortex 20-30min, spin, sonicate 5 min. Repeat steps 1 and 2.
16. Combine and vortex the extracted digests, spin and Speed Vac to reduce volume to 10  $\mu\text{L}$ .
17. Either proceed with C18 ZipTip (Millipore) cleanup or analyze with LC-MS. Add 2-5  $\mu\text{L}$  of 5% formic acid. When analyzing low levels of protein, concentrate the peptides by eluting from ZipTips using 3 $\mu\text{L}$  of elution solution into a clean 0.65 mL siliconized tube.

Adapted from UCSF MS Facility: <http://ms-facility.ucsf.edu/ingel.html>